



# DIAGNOSTIC INSTRUMENTS



**01-101-00**  
Large Face Aneroid Model  
Blood pressure unit



**01-102-00**  
Mercury Sphyg. Desk Model  
Stethoscope Blood pressure unit



**01-103-00**  
Standard Aneroid Adult Cuff  
Stethoscope Blood pressure unit



**01-104-00**  
Standard Aneroid Child Cuff  
Stethoscope Blood pressure unit



**01-105-00**  
Standard Aneroid Large Adult  
Cuff Blood pressure unit



**01-106-00**  
Fanny Pack Aneroid Latex Free with  
instruments Blood pressure unit



**01-107-00**  
Fanny Pack Aneroid Latex Free  
without instruments Blood pressure unit



**01-108-00**  
Home Blood Pressure Kit  
Blood pressure unit



**01-109-00**  
Blood Pressure Unit for regular  
and large arm



**01-110-00**  
Wrist blood  
pressure monitor



Dual Head Stethoscope  
**01-111-00** Black tubing  
**01-112-00** Red tubing  
**01-113-00** Burgandy tubing



**01-114-00**  
Dual Head Stethoscope  
Red tubing



**01-115-00**  
Dual Head Stethoscope  
Burgandy tubing



Single Head Stethoscope  
**01-116-00** Black tubing  
**01-117-00** Red tubing  
**01-118-00** Burgandy tubing



**01-119-00**  
Single Head Stethoscope  
red tubing



**01-120-00**  
Dual Head Stethoscope  
Burgandy tubing



**01-121-00**  
Stainless Steel  
Stethoscope



**DUPLEX**  
Stethoscope  
01-122-00  
super lightweight



**DUPLEX-BABY**  
Stethoscope  
01-123-00  
super lightweight



**ANAESTOPHON**  
Stethoscope  
01-124-00  
super lightweight



Clar Mead Mirror  
komplett  
6 Volt  
01-125-00



Ziegler  
01-126-00  
Ø 90mm  
with head band without  
protection cap



Ziegler  
01-127-00  
Ø 90mm  
with head band with  
protection cap



01-128-00  
Ø 50mm



**PINARD WOODEN**  
Stethoscope  
01-129-00  
18 cm




**PINARD ALUMINIUM**  
Stethoscope  
01-130-00  
15 cm



**WEBER**  
Aesthesiometer  
01-131-00  
9 cm

**Blood Lancets**




**FRANCKE**  
01-132-00

**FRANCKE**  
01-133-00

**CHEST PIECES**



**BOWLES**  
01-134-00



**FORD**  
01-135-00



**LIGHT WEIGHT**  
01-136-00



**SPRAGUE-BOWLES**  
01-137-00



**BOWLES**  
01-138-00



**FORD**  
01-139-00



**BELL-END**  
01-140-00



**SPARE RUBBER**  
01-141-00



**DIAPHRAM**  
01-142-00



**SCHIMMELBUSCH**  
01-143-00 ADULT  
01-144-00 CHILD

**INHALERS**

**AIRWAYS**



**YANKAUER**  
01-145-00 ADULT  
01-146-00 CHILD



**CONNELL**  
01-147-00 ADULT  
01-148-00 CHILD



**MAGILL**  
Catheter Introducing  
Forcep  
01-149-17 17cm



**MAGILL**  
Catheter Introducing  
Forcep  
01-150-20 20cm



**MAGILL**  
Catheter Introducing  
Forcep  
01-151-25 25cm

## PERCUSSION HAMMERS

Perkussionshämmer | Martillos de percusión  
Marteaux-percuteurs | Martelli percussori



**NEURO-DUOFLEX Hammer**  
01-152-22 22cm



**DEJERINE Hammer**  
01-153-20 20cm



**DEJERINE Hammer**  
01-154-20 20cm



**TROEMNER Hammer**  
01-155-14 145 g  
01-156-17 175 g  
01-157-23 235 g  
24.5cm  
01-158-00 spare rubber



**TROEMNER-SIMPLEX Hammer**  
01-159-21 21cm



**BUCK Hammer**  
01-160-18 18.5cm  
01-161-20 20cm



**BUCK HAMMER**  
01-162-18 18cm  
01-163-20 20cm



**TAYLOR Hammer**  
01-164-16 16cm  
01-165-18 18cm  
01-166-20 20cm



**TAYLOR Hammer**  
01-167-18 18cm



**Babinsky Hammer**  
01-168-22 22cm



**Babinsky Hammer**  
01-169-22 22cm




**Rabiner Hammer**  
01-170-25 25.5cm  
with needle & brush



**National Hospital type  
Percussion Hammers**  
Length diameter  
01-171-50 350mm, 50mm  
01-172-60 350mm, 60mm  
01-173-70 350mm, 70mm



**Berliner Hammer**  
01-174-20 20cm




**Traube Hammer**  
01-175-16 16cm



**Wartenberg Pinwheel**  
01-176-19 19cm

## Conventional Laryngoscope Set

### MC-INTOSH LARYNGOSCOPE



#### BLADE

Size	A	B	C	Working Length
0	80mm	58mm	30mm	55mm
1	97mm	76mm	36mm	70mm
2	111mm	88mm	36mm	78mm
3	133mm	112mm	36mm	108mm
4	154mm	131mm	38mm	133mm

#### Mc-Intosh Laryngoscope

Plastic case, Standard illumination. Manufactured from AISI 304 Stainless Steel Matt Type, Standard Illumination 2.5 V.

#### Blades Availability:

- ★ Matt Finish
- ★ Mirror Finish
- ★ Sand Blast Finish
- ★ Pls, confirm in your order list

#### Mc-INTOSH LARYNGOSCOPE SET STANDARD ILLUMINATION

##### 01-177-00

##### PEDIATRIC SET

Set of 3 Blades

Fig. 0,1,2

Supplied in plastic case with 2 spare blade

##### 01-178-00

##### STANDARD SET

Set of 3 Blades

Fig. 1,2,3

Supplied in plastic case with 2 spare blade

##### 01-179-00

##### ADULT SET

Set of 4 Blades

Fig. 1,2,3,4

Supplied in plastic case with 2 spare blade

##### 01-180-00

##### ADULT SET

Set of 5 Blades

Fig. 0,1,2,3,4

Supplied in plastic case with 2 spare blade

Fig. 0

01-181-00

Fig. 1

01-182-01

Fig. 2

01-183-02

Fig. 3

01-184-03

Fig. 4

01-185-04



01-186-20  
Spare Bulb



01-187-10  
Laryngoscope Battery  
Handle for Adult  
C type 2 Batteries




Mc-Intosh Laryngoscope Set  
Standard Illumination  
01-188-00



## Conventional Laryngoscope Set

### MILLER LARYNGOSCOPE



**BLADE**

Size	A	B	C	Working Length
0	75mm	53mm	26mm	55mm
1	111mm	80mm	26mm	80mm
2	151mm	134mm	28mm	130mm
3	188mm	167mm	28mm	170mm
4	203mm	180mm	28mm	180mm

#### Mc-Intosh Laryngoscope

Plastic case, Standard illumination.

Manufactured from AISI 304 Stainless Steel.

#### Blades Availability:

- \* Matt Finish
- \* Mirror Finish
- \* Sand Blast Finish
- \* Pls, confirm in your order list

### MILLER LARYNGOSCOPE SET OF CONVENTIONAL STANDARD ILLUMINATION

#### 01-189-00 PEDIATRIC SET

Set of 3 Blades  
Fig. 0,1,2  
Supplied in plastic case with 2 spare blade

#### 01-190-00 STANDARD SET

Set of 3 Blades  
Fig. 1,2,3  
Supplied in plastic case with 2 spare blade

#### 01-191-00 ADULT SET

Set of 4 Blades  
Fig. 1,2,3,4  
Supplied in plastic case with 2 spare blade

#### 01-192-00 ADULT SET

Set of 5 Blades  
Fig. 0, 1,2,3,4  
Supplied in plastic case with 2 spare blade

Fig. 0  
01-193-00

Fig. 1  
01-194-01

Fig. 2  
01-195-02

Fig. 3  
01-196-03


Fig. 4  
01-197-04

01-198-20  
Spare Bulb

01-199-10  
Laryngoscope Battery  
Handle for Adult C type  
2 batteries

Laryngoscope Set of conventional Miller  
Standard Illumination  
01-200-00

## MC-INTOSH FIBER OPTIC LARYNGOSCOPE BLADE



BLADE

Size	A	B	C	Working Length
0	77mm	53mm	14mm	55mm
1	95mm	72mm	43mm	70mm
2	111mm	86mm	44mm	78mm
3	133mm	109mm	44mm	108mm
4	151mm	129mm	44mm	133mm

### Mc-Intosh Laryngoscope Fiber Optic

Halogen illumination. Manufactured from AISI 304 Stainless-Steel. Halogen Illumination 2.5 V. The Bulb is inside the handle in order to simplify the sterilization of the blades.

#### Blades Availability:

- \* Matt Finish
- \* Mirror Finish
- \* Sand Blast Finish
- \* Pls, confirm in your order list

## Mc-INTOSH FIBER OPTIC LARYNGOSCOPE SET HALOGEN ILLUMINATION

**01-201-00 PEDIATRIC SET**  
Set of 3 Blades  
Fig. 0,1,2  
Supplied in plastic case with 2 spare blade

**01-202-00 STANDARD SET**  
Set of 3 Blades  
Fig. 1,2,3  
Supplied in plastic case with 2 spare blade

**01-203-00 ADULT SET**  
Set of 4 Blades  
Fig. 1,2,3,4  
Supplied in plastic case with 2 spare blade

**01-204-00 ADULT SET**  
Set of 5 Blades  
Fig. 0,1,2,3,4  
Supplied in plastic case with 2 spare blade



Fig. 0  
01-205-00



Fig. 1  
01-206-01



Fig. 2  
01-207-02



Fig. 3  
01-208-03



Fig. 4  
01-209-04



01-210-20  
Laryngoscope  
Bulb Small



01-211-22  
Laryngoscope  
Bulb Large




01-212-10  
Laryngoscope Battery  
Handle for Adult  
C type 2 batteries



01-213-00  
Mc-Intosh Fiber Optic Laryngoscope Set  
Halogen Illumination

## MILLER FIBER OPTIC LARYNGOSCOPE



### BLADE

Size	A	B	C	Working Length
0	75mm	53mm	26mm	55mm
1	111mm	80mm	26mm	80mm
2	151mm	134mm	28mm	130mm
3	188mm	167mm	28mm	170mm
4	203mm	180mm	28mm	180mm

### MILLER FIBER OPTIC LARYNGOSCOPE

Halogen illumination. Manufactured from AISI 304 Stainless-Steel. Halogen Illumination 2.5 V. The Bulb is inside the handle in order to simplify the sterilization of the blades.

#### Blades Availability:

- ★ Matt Finish
- ★ Mirror Finish
- ★ Sand Blast Finish
- ★ Pls, confirm in your order list

### MILLER FIBER OPTIC LARYNGOSCOPE HALOGEN ILLUMINATION

**01-214-00 PEDIATRIC SET**  
Set of 3 Blades  
Fig. 0,1,2  
Supplied in plastic case with 2 spare blade

**01-215-00 STANDARD SET**  
Set of 3 Blades  
Fig. 1,2,3  
Supplied in plastic case with 2 spare blade

**01-216-00 ADULT SET**  
Set of 4 Blades  
Fig. 1,2,3,4  
Supplied in plastic case with 2 spare blade

**01-217-00 ADULT SET**  
Set of 5 Blades  
Fig. 0,1,2,3,4  
Supplied in plastic case with 2 spare blade



Fig. 0  
01-218-00



Fig. 1  
01-219-01



Fig. 2  
01-220-02



Fig. 3  
01-221-03



Fig. 4  
01-222-04



**01-223-20**  
Laryngoscope Bulb Small



**01-224-22**  
Laryngoscope Bulb Large



**01-225-10**  
Laryngoscope Battery Handle for Adult  
C type 2 Batteries



**01-226-00**  
MILLER FIBER OPTIC LARYNGOSCOPE

## Laryngoscope Flexible (Mc-Coy)

### CONVENTIONAL SET

**01-227-00 Flexible Laryngoscope Conventional Set of 3 blades**  
Fig. 2,3,4

**01-228-00 PEDIATRIC SET**  
Handle With Blade Fig. 2  
Supplied in plastic case

**01-229-00 STANDARD SET**  
Handle With Blade Fig. 3  
Supplied in plastic case

**01-230-00 ADULT SET**  
Handle With Blade Fig. 4  
Supplied in plastic case  
Halogen Illumination  
Manufactured From  
Stainless Steel High  
grade 304 material.



**01-231-10 'C' type battery handle**  
Blades Size:  
**01-232-02** Fig. 2  
**01-233-03** Fig. 3  
**01-234-04** Fig. 4

KUMMAS  
CORPORATION

**01-235-00 Flexible Laryngoscope Conventional Set of 3 blades**  
Fig. 2,3,4

**01-236-00 PEDIATRIC SET**  
Handle With Blade Fig. 2  
Supplied in plastic case

**01-237-00 STANDARD SET**  
Handle With Blade Fig. 3  
Supplied in plastic case

**01-238-00 ADULT SET**  
Handle With Blade Fig. 4  
Supplied in plastic case

Halogen Illumination  
Manufactured From  
Stainless Steel High  
grade 304 material  
packed in green plastic  
box.



**01-243-00**



**01-239-10 'AA' type battery handle**  
Blades Size:  
**01-240-02** Fig. 2  
**01-241-03** Fig. 3  
**01-242-04** Fig. 4



Two laryngeal mirrors: n. 4 and 5



Bent-arm illuminator: pin contact fitting.



Specula



With 2 spare Bulbs

### Professional Otoscope & Ophthalmology Set ENT Medical Diagnostic Set

01-244-00

2.5 V Standard illumination for use with 2 batteries C type.

Including:

Otoscope Head: Luer Lock fitting.

Chrome plated.

Fix lens 4 x in plastic mount can be removed for introduction of instruments.

Specula sizes: 2.5, 3.5 & 4.5 mm Included,

Ophthalmoscope head: Luer Lock fitting.

Lens selection wheel permits choice of corrective lenses from 0 to - 20 and 0 to +20 diopters.

Single large circle aperture.

Bent-arm illuminator: pin contact fitting.

Nasal speculum: Chrome plated.

Two laryngeal mirrors: n. 4 and 5.

Battery handle: 28mm. For 2 cell C-batteries.

Chrome plated. For adjustment of all illumination.

Replacement bulb: 2.5 V

Diagnostic set is supplied in Quality Plastic Case.



Otoscope Head



Nasal speculum: Chrome plated.



Ophthalmoscope head



ENT Specula



Two laryngeal mirrors: n. 4 and 5



Bent-arm illuminator: pin contact fitting.

### Universal Otoscope & Ophthalmology Set ENT Medical Diagnostic Set

01-245-00

2.5 V Standard illumination for use with 2 batteries type C.

Including:

Otoscope Head: Luer Lock fitting.  
Chrome plated.

Fix lens 4 x in plastic mount can be removed for introduction of instruments.

Specula sizes: 2.5, 3.5 & 4.5 mm Included,  
Ophthalmoscope head: pin contact fitting.

Lens selection wheel permits choice of corrective lenses from 0 to - 20 and 0 to +20 diopters.

Single large circle aperture.

Bent-arm illuminator: Luer Lock fitting.

Nasal speculum: Chrome plated.

Two laryngeal mirrors: n. 4 and 5.

Battery handle: 28mm. For 2 cell C-batteries.

Chrome plated. For adjustment of all illumination.

Replacement bulb: 2.5 V

Diagnostic set is supplied in Quality Plastic Case.



Nasal speculum: Chrome plated.



ENT Specula



Otoscope Head



Ophthalmoscope Head



Handle



### DIAGNOSTIC PENTOSCOPE SET 01-246-00

2.5 V Standard illumination for use with 2 batteries type AA.

Consisting:

Universal Handle thread type fit for all Pentoscope Heads like.

- Otoscope Head.
- Ophthalmoscope Head.
- Throat Branch Head.
- Mouth Mirrors Head.

It have two exploring mirrors No. 4 & 5 mm.

Its spare 06 disposable Plastic ear Canulas are with sizes 2 x 2.50 mm, 2 x 3.50 mm & 2 x 4.5 mm. On special order it can be as 2 x 1.5 mm for new born babies 2 x 2.5 mm for pediatric and 2 x 4 mm for adults, so it means total 06 canula also available in a plastic tube packing.

Sturdy adjustment mount bracket

2 spare 2.5 V bulb.

Pentoscope set battery handle is operating with two AA size batteries.

This Pentoscope set Packed in our new beautiful Blue Colour plastic case in same blue inner packing material. Spare parts also available of any item.



With 2 spare Bulbs



Mouth Mirrors Head



Standard illumination



# KUMMAS

C O R P O R A T I O N



Otoscope Head



Sturdy Adjustment Mount Bracket



Plastic ear Canulas



**Standard Otoscope**  
01-247-00

This otoscope is 2.5 V standard illumination, handle is powder coated as well otoscope head too. Handle and head is Pin Screw lock fitting. Head of otoscope is with 6x magnifying glass. 03 Reusable speculas with sizes 2.5, 3.5 & 4.5 mm. Battery handle for 2 'c' batteries. Rheostat for adjustment of illumination. Supply in Nice Wooden Box Packing



**Standard Otoscope**  
01-249-00

This otoscope is 2.5 V standard illumination, handle is chrome plated as well otoscope head too. Handle and head is Luer / Faster lock fitting. Head of otoscope is with 6x magnifying glass. 03 Reusable speculas with sizes 2.5, 3.5 & 4.5 mm. Battery handle for 2 'c' batteries. Rheostat for adjustment of illumination. Supply in Black Plastic Case Packing



**Standard Otoscope**  
01-248-00

This otoscope is 2.5 V standard illumination, handle is chrome plated as well otoscope head too. Handle is Pin Screw lock fitting. Head of otoscope is with 6x magnifying glass. 03 Reusable speculas with sizes 2.5, 3.5 & 4.5 mm. Battery handle for 2 'c' batteries. Rheostat for adjustment of illumination. Supply in Black Plastic Case Packing



**Mini Otoscope Fiber Optic**  
01-250-00

2.5 V fiber optic halogen light. Otoscope head: pin contact fitting. Fixed lens 6X in plastic mount may be removed for introduction of instruments. Reusable specula 2, 3 and 4 mm. Included Battery Handle Made of Brass, Ionized Plated handle with Push button for 2 AA Batteries. Replacement bulb: 2.5 v the diagnostic set is supplied in Black Pouch. Available in multi colour.



**Mini Otoscope Fiber Optic**  
01-251-00

2.5 V fiber optic halogen light. Otoscope head: pin contact fitting. Fixed lens 6X in plastic mount may be removed for introduction of instruments. Reusable specula 2, 3 and 4 mm. Included Battery Handle Made of Brass, Ionized Plated handle with Push button for 2 AA Batteries. Replacement bulb: 2.5 v the diagnostic set is supplied in Black Pouch. Available in multi colour.







**Standard Otoscope + Ophthalmoscope Set**  
**01-252-00**

This Set is with 2.5 V standard illumination, handle is chrome plated as well otoscope head too. Handle and head is Luir / Faster lock fitting. Head of otoscope is with 6x magnifying glass. 03 Reusable speculas with sizes 2.5, 3.5 & 4.5 mm. A Full Functional Ophthalmoscope compact durable, light weight. Standard bulb and precision system 4 different apertures selection range of lenses -20 to +20 bright white light and brilliant image A Full Functional Ophthalmoscope Battery handle for 2 'c' batteries. Rheostat for adjustment of illumination. Supply in Black Plastic Box Packing.



Otoscope Head



ENT Specula



Ophthalmoscope head



ENT Specula



'C' type Battery handle

C O R P O R A T I O N



**Standard Ophthalmoscope Set**  
**01-253-00**

A Full Functional Ophthalmoscope compact durable, light weight. Standard bulb and precision system 4 different apertures selection range of lenses -20 to +20 bright white light and brilliant image Packed in Black Plastic Case.



Ophthalmoscope head



'C' type Battery Handle



**Standard Otoscope Set  
01-254-00**

this otoscope is with 2.5 V standard illumination, handle is chrome plated as well otoscope head too. Handle and head is Luer / Faster lock fitting. Head of otoscope is with 6x magnifying glass. 03 Reusable speculas with sizes 2.5, 3.5 & 4.5 mm. Battery handle for 2 'c' batteries. Rheostat for adjustment of illumination. Supply in Black Pouch Packing



**Standard Otoscope Set  
01-255-00**

this otoscope is with 2.5 V standard illumination, handle is chrome plated as well otoscope head too. Handle and head is Luer / Faster lock fitting. Head of otoscope is with 6x magnifying glass. 03 Reusable speculas with sizes 2.5, 3.5 & 4.5 mm. Battery handle for 2 'c' batteries. Rheostat for adjustment of illumination. Supply in Black Pouch Packing




**Otoscope + Ophthalmoscope Set  
01-256-00**

This otoscope is with 2.5 V standard illumination, handle is chrome plated as well otoscope head too. Handle and head is Luer / Faster lock fitting. Head of otoscope is with 6x magnifying glass. 03 Raisable speculas with sizes 2.5, 3.5 & 4.5 mm. A Full Functional Ophthalmoscope compact durable, light weight. Standard bulb and precision system 4 different apertures selection range of lenses -20 to +20 bright white light and brilliant image A Full Functional Ophthalmoscope Battery handle for 2 'c' batteries. Rheostat for adjustment of illumination. Supply in Black Pouch Packing.



**Mini Ophthalmoscope  
01-257-00**

A Full Functional Ophthalmoscope.  
compact durable, light weight.  
Standard bulb and precision system  
4 different apertures selection  
range of lenses -20 to +20  
bright white light and brilliant image  
Made of plastic Material.



**Mini Ophthalmoscope  
01-258-00**

A Full Functional Ophthalmoscope.  
compact durable, light weight.  
Standard bulb and precision system  
4 different apertures selection  
range of lenses -20 to +20  
bright white light and brilliant image  
made of plastic Material.



**Standard Ophthalmoscope  
01-259-00**

A Full Functional Ophthalmoscope.  
compact durable, Standard bulb and  
precision system  
4 different apertures selection  
range of lenses -20 to +20  
bright white light and brilliant image  
Handle made of Stainless Steel Material.



**GRADUATED TAPE MEASURER**  
**01-260-02**  
2 m long (1.5 m) graduated in inch & cm markings



**01-261-15** 15cm  
**01-262-30** 30cm  
cm & mm markings



**JAMESON CALIPER**  
**01-263-08**  
8cm



**TOWNLEY CALIPER**  
**01-264-10**  
10cm

**KUMMAS**  
CORPORATION



**VERNIER CALIPER**  
**01-265-15**  
15cm

**VERNIER CALIPER**  
**01-266-15**  
15cm



**MÖELTGEN Flexometer**  
**01-267-24**  
24cm



**CARROLL POCKET GONIOMETER**  
**01-268-14**  
14.5cm



**BERTILLON MEASURING DEVICE**  
**01-269-32**  
32cm



**CASTROVIEJO-EPKER  
Caliper**  
01-270-09 9 cm Str.  
Measuring range, 40 mm  
01-271-09 9 cm Cvd.  
Measuring range, 40 mm



**CASTROVIEJO-EPKER  
Caliper**  
01-272-18 18 cm  
Measuring range  
40 mm, straight



**CASTROVIEJO-EPKER  
Caliper**  
01-273-18 18 cm Angled  
Measuring range  
40 mm



**CASTROVIEJO  
01-274-18 18cm  
Bone Caliper  
Straight**

KUMMAS CORPORATION



**CASTROVIEJO  
01-275-18 18cm  
Bone Caliper  
Curved**



**CASTROVIEJO  
01-276-08  
8 cm  
Eye View Caliper  
Measuring range  
20 mm**



**Tonometer  
01-277-00  
Set**

